




For Lifting & Tightening


Forehead Wrinkle




Crow's Feet




Marionette Lines




Fine Wrinkles




Lip Lines




Cheek Definition




Neck Folds




Collagen synthesis declines in the dermis over time.



Focused ultrasound shrinks skin by delivering energy directly into it.




Stimulated by focused energy, skin collagen regenerates.




By regenerating collagen, the skin is firmed and tightened, which eliminates wrinkles.


For Slimming & Contouring




Fat cells and cellulite is formed within the body.



The smasfocus X energy destroys the targeted fat cells.



As part of the body’s natural healing response, cells regenerate and join together to fill the gap.



The fat layer is naturally eliminated from the body, resulting in smooth and well-contoured body lines.

Thorax

Axilla


Arms

Flanks

Abdomen

Thighs

Single-line Handle



NORMAL MODE

MICRO PULSE MODE



Change the Mode with One Click

7 Cartridges to Meet Various Body Part Needs

10000 SHOTS

FACE



1.5mm, 2.0mm, 3.0mm, 4.5mm

BODY



6.0mm, 9.0mm, 13.0mm

Multi-line Handle



Reduced Treatment Time

The groundbreaking S-shaped working mode can reduce nearly twice the treatment time.

Various Cartridges Up to 8 Types

20000 SHOTS

FACE



1.5mm, 3.0mm, 4.5mm

BODY



6.0mm, 8.0mm, 10.0mm, 13.0mm, 16.0mm

Booster Handle

Frequency: 1-8hz

(90°–360° Adjustable)

NORMAL

CIRCULAR

MICRO CIRCULAR

Double lifting effects with vmax overlapping treatment after the linear procedure.

Various Cartridges Up to 3 Types

60000 SHOTS

FACE



1.5mm,

3.0mm,

4.5mm

Making the handpiece smaller in size

Generating energy in a circular motion



Slim and light cartridges with adjustable penetration depths. Convenient and safe treatment for curved and narrow areas with better visibility on the treatment spot.

Private Handle

FOR VAGINAL TIGHTENING


Laser vs HIFU

- Collagen Rejuvenation
- Tightening


More Stronger

More Deeper


More Safer



10000 SHOTS




3.0mm, 4.5mm



AUTOMATED OPERATION

5–90 Rotation Angle Adjustable





15.6-INCH HD SCREEN

Ultra-thin border ultra-large display screen ratio



Aviation-grade Power Supply

Built-in instant response sensor significantly enhances safety performance.



01 Mobile Response Sensor

Once the piezoelectric ceramic stops, handle immediately shuts down to prevent harm.

02 Allegro MicroSystems Stepper Motor

American imported motor with extended service life.

03 Japanese Piezoelectric Ceramics

Ensures continuous and stable energy output.

Before & After



Specifications

Power	1200 W
Power Frequency	4mhz–7mhz
Screen	15.6 Inches
Width	1–2.5mm Adjustable
Length	8–25mm Adjustable
Dot Pitch	1.0–2.5mm Adjustable
Energy Output	0.1–2J/cm2